RITA radioactivity TLC analyser





General description

RITA – the radioactivity thin-layer-chromatography detector – is using a linear analyser detector, which is sensitive over the entire chromatogram trace from the start to the front. The entrence window can be open for 3H-applications or closed for 14C or any other beta[–] radiation emitting nuclides. Many gamma radiation emitting nuclides are emitting beta radiation as well. The pure gamma emitting nuclides can be detected as well, because they produce some Compton electrons.

RITA has 2 simultaneous sensing electrodes: the delay line determines the location of the counted event, the pulse height of the counting wire is used for the electronic collimation. Thus high energy beta radiation can be detected with high resolution using a narrow collimator window. RITA can measure many single traces one after the other. The Rita detector is elevated after the end of the measurement and moved automatically to the next trace position. Multi traces can be dispayed in 3-dimensional presentation or a 2-dimensional distribution over the entire 200 x 200 mm TLC plate can be calculated and presented.

Counting gas is a non-inflammable mixture of 90% Argon and 10% Methane, which is traded as P10 gas. The gas consumption is very low.

2 TLC plates of 200 x 200 mm each can be placed in Rita at the same time.

A metal/glass cover closes Rita completely and protects the environment as well as Rita itself from electronic frequency emissions (CE).

Applications

- β- radioactivity detector for Thin-Layer-Chromatography
 (β+)
 - nuclides 3H, 14 C etc nuclides 18F, 11 C etc

external standardisation

ultra high sensitivity

Features

open window detection for 3H

antistatic protection grid for open window operation

closed window detection for 14C etc

very high sensitivity

live display of measurement

peak integration, backgound subtraction

manual or automatic

limit of detection calculation

Ordering information

- **02110001** RITA * B-TLC anlyser for single and multiple traces 2 and 3 dimensional presentation
- 02000105 RITA* homogeneity +resolution standard 14C source 22140000dpm 200 x 10 mm with holder
- **02110010** RITA*IQ/OQ 14C standard for check of homogeneity + resolution 18 spots homogeneity test 4 spots sensitivity test 1-3-mm resolution test

counting gas and valve are expected on-site

& raytest



400 5mm 4 a 3, ... resolution test



uniformity test

Те	ch	ni	cal	d	ata
	~		- ui		utu

detector:	gas flow proportional counter		
counting gas:	P10 90% Argon, 10% Methane		
gas consumption:	11/min		
detection area: single step	200 mm long, 20 mm wide		
detection window:	open or closed, 0,9 mg/cm2		
diaphragm:	3 - 20 mm wide magnetic attached		
resolution: 3H	< 1 mm < 2 mm		
sensitivity: 3H 14C	1000 dpm in 10 min 100 dpm in 10 min		
background:	80 dpm/200 mm		
energy discrimination:	electronic collimator		
traces:	1-80		
display	live single chromatogram		
evaluation:	manual or automatic peak integration background subtraction limit of detection		
dimension (wxhxd):	700 x 320 x 560 mm		
weight:	45 kg		
power supply:	110-230 V, 20 VA		
operating conditions:	10-40°C max 70% r.H.		

Sole Distributor in China and Hong Kong : NIKYANG Enterprise Limited

The second secon Shanghai • 86-21 6351 1828

info@nikyang.com

Guangzhou • 86-20 8329 2451 Beijing • 86-10 6527 8522

www.nikyang.com

