

SOTAX MPS Media Preparation Station



- Efficient preparation of dissolution media within a mobile space-saving unit
- Open system, easy to clean, observe and validate
- Exact preheating of the media to test temperature and reliable degassing
- Continuous checks of selected temperature
- Fast and accurate volumetric filling of test vessels at 2000 ml/min
- Use of concentrates with on-line preparation or by sequential addition
- Option of gravimetric self-calibration with built-in balance and automatic report generation



Media Preparation Station for the efficient preparation of dissolution media (degassing, heating and dosing of media)





Self-calibration

Technical data

Preparation batches 10 or 20 litres.

switchable

30-40 °C Temperature range Preparation time 10 / 20 min Medium feed rate 2000 ml/min

Dosing accuracy 50-80 ml +/-0.8 ml,

80-1000 ml +/-1% 230 V, 50 Hz or

Voltage 115 V, 60 Hz

Power consumption 2000 VA 0.2 bar Vacuum

Connecting hoses for:

Medium supply 8/10 dia. Cleaning water 3/4 " Waste 8/10 dia.

Media usual dissolution

media (no organic

solvent)

Width, depth, height 450, 690, 930 mm

Order information

Part # Description

SOTAX MPS Media 8200

Preparation Station

8346 Self-calibration

option

Key hardware features

When USP 1,2,4,5 & 6 dissolution tests are carried out, the preparation of the media and its degassing, tempering and accurate metering into the test vessels involves a considerable amount of work. The use of an automatic media preparation is one step in the automation of dissolution tests in the laboratory.

The SOTAX MPS Media Preparation Station permits the efficient preheating and degassing of the dissolution media as well as fast and accurate filling. In laboratories where many short tests are carried out or where several dissolution units are used, these can be filled with media quickly and easily. With the use of concentrates there is the possibility to gain even more efficiency.

Description

The SOTAX MPS Media Preparation Station is very compact and can be easily moved around the laboratory. A volume of media of 10 or 20 litres may be prepared which is sufficient for up to three dissolution tests. The unit was designed to be easy to use, clean, program and validate.

Mode of operation

The SOTAX Media Preparation Station has a Plexiglas media container and works with a liquid volume of up to 20 litres. The container is under vacuum. During this process the media is degassed and simultaneously heated with the SOTAX Thin Film Technology™. Calibration to the desired temperature is effected with a second heater placed inside the container.

During the preparation process the degassing and temperature-adjusting functions are automatically controlled.

The volume of media can be selected via the front key pad and is metered from the media container into the test vessel using a highly accurate and fast delivery pump with 2000 ml/min flow rate on a single channel. The dissolution tester is filled with media and ready to use in approximately three minutes. The accuracy is >99.5 % fulfilling the USP requirements. Before a different medium is prepared, an automatically controlled emptying and cleaning function is set in motion. The glass vessel and the entire tubing and valve system are thoroughly flushed to avoid contaminating the subsequent media. The system's open architecture allows for easy visible inspection of the cleaning results.

The SOTAX Media Preparation Station is equipped for sequential metering of deaerated deionised water and media concentrates. The dilution ratio may be programmed via the front key pad. Using the SOTAX Media Preparation Station in this mode eliminates cleaning of the station between media

A sophisticated electronic system automatically controls and monitors all functions during the preparation and cleaning process. Key data is shown on the display and may be printed by connecting a printer to the serial port.

Self-calibration (option)

Self-calibration system with built-in balance for ease of validation. The menu-guided process requires less than one minute and the result may be printed as validation documentation. Up to 5 validation readings are possible.

Sole Distributor in China and Hong Kong: NIKYANG Enterprise Limited

Hong Kong • 852-2569 2154

Guangzhou • 86-20 8329 2451

Shanghai • 86-21 6351 1828

• 86-10 6527 8522 Beijing





